

ESC Endoscope Storage Cabinet

Troubleshooting Alarms

- [Log In Failed](#) - Re-enter or scan USER ID. If this fails contact the administrator to either amend your pin or USER ID.
- [Endoscope Not Found](#) - This is indication that the ENDOSCOPE ID is blank and holds no relevant scope data. This can only be assigned by an **administrator**.
- [Airflow Alarm](#) - Indication that the cabinet filters are blocked or require replacement.
- [Close The Doors](#) - The cabinet will not proceed with any further instructions or monitoring until the doors are closed.
- [Air Pump Failure 1](#) - Contact Medivators Technical Support. Do not use the left side of the cabinet, and reprocess the scopes stored here.
- [Air Pump Failure 2](#) - Contact Medivators Technical Support. Do not use the right side of the cabinet, and reprocess the scopes stored here.
- [Air Pump Failure 1 & 2](#) - Both pumps have failed. Contact Medivators Technical Support. Do not use the the cabinet and reprocess the scopes stored here.
- [Service Due](#) - displayed 180 days after the unit has been serviced to warn that a service is due in 28 days. Failure to have the cabinet serviced regularly may result in storage validity not being verified and the scopes may be unsafe to use directly from cabinet.