

ESC Endoscope Storage Cabinet

Daily Use Protocol

Loading an Endoscope

- NOTE: CHOOSE AN APPROPRIATE METHOD OF HAND DECONTAMINATION ACCORDING TO YOUR LOCAL POLICY BEFORE HANDLING ENDOSCOPES.
- Once the endoscope has been reprocessed it should be placed into the cabinet without delay.
- Scan USER ID or log in manually. If screen is dark, touch screen to animate it.
- Select ENDOSCOPE IN on screen
- Choose hook number
- Scan ENDOSCOPE ID or select OK to unlock door.
- Place scope onto chosen hook.
- Once the scope is in place, ensure all channels are appropriately connected and the channel separator is in place.
- Connect the DSD-hookup to the cabinet by inserting the male CPC connector into the cabinet manifold until a “click” is heard.
- Close door and press OK, a print out of the scope details will then be produced.
- Once in the cabinet the scope will remain viable for up to the preset storage time of depending on facility hang-time policies.

Removing an Endoscope

- Scan USER ID or or log in manually.
- Select ENDOSCOPE OUT on screen
- Choose which door to open.
- Disconnect the DSD hookup from the cabinet and remove endoscope.
- Scan the ENDOSCOPE ID, close the door and press OK to confirm what scope has been removed.
- A printout will then be produced with scope status, storage details, and user load/unload information.