CAT. No. ML02-0133

Contenu Net:
3,8 L (1 gal US)
(3,8 litres)

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Low foaming, neutral pH, bactericidal detergent specifically formulated for use with medical instruments

Fabriqué aux É.U. pour:

MEDIVATORS Inc.
14605 28th Avenue North
Minneapolis, MN 55447 U.S.A.
Toll Free: 1.800.444.4729

www.medivators.com

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VANTAGE ENDOSCOPY

90 Gough Road, Unit 1
Markham, ON L3R 5V5O
1-866-677-4121

CAT. No. ML02-0133

Net Content:
1 U.S. Gallon
(3.8 Liters)

KEEP OUT OF REACH OF CHILDREN

CAUTION

DIRECTIONS FOR USE: For cleaning of fully immersible endoscopes, related accessories, surgical instruments, and other apparatus where blood, mucus, protein or other hard to remove soils are encountered, use INTERCEPT at 1/3 oz./gal of water (0.25% use concentration) **with one full stroke of the hand-pump (1 oz.) to three gallons of water.** If visible organic matter is present, use INTERCEPT at 2/3 oz./gal of water (0.5% concentration) **which is two full strokes of the hand-pump (2 oz.) to three gallons of water.** For best manual cleaning results, mix INTERCEPT with cool to warm water (20°C-35°C) (38°F – 95°F) and ensure a minimum contact time of (1) minute. Rinse all surfaces and internal channels thoroughly with water. For use in automated washers, follow washer manufacturer's recommendations. Do not mix INTERCEPT with bleach.

PRECAUTIONARY STATEMENTS: Use of personal protective equipment including appropriate clothing, gloves, and safety glasses must be worn. Refer to the material Safety Data Sheet (MSDS) for additional information.

STATEMENT OF PRACTICAL TREATMENT: In case of contact with eyes, flush with copious amounts of water for 15 minutes. In case of ingestion, DO NOT induce vomiting. Contact CHEMTREC for emergency medical treatment at 1-800-424-9300.

STORAGE & DISPOSAL: Store at room temperature 15°C-25°C (59°F-77°F). INTERCEPT has a two (2) year shelf life from date of manufacture. Discard unused solution according to local, state and federal regulations. DO NOT reuse container.

Read MSDS for additional precautionary statements

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© 2013 MEDIVATORS Inc. 67199-666/C