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**VANTAGE ENDOSCOPY**90 Gough Road, Unit 1  
Markham, ON L3R 5V5O  
1-866-677-4121**CAT. No. ML02-0133****Net Content:  
1 U.S. Gallon  
(3.8 Liters)****KEEP OUT OF REACH OF CHILDREN****CAUTION**

**DIRECTIONS FOR USE:** For cleaning of fully immersible endoscopes, related accessories, surgical instruments, and other apparatus where blood, mucus, protein or other hard to remove soils are encountered, use INTERCEPT at 1/3 oz./gal of water (0.25% use concentration) **with one full stroke of the hand-pump (1 oz.) to three gallons of water.** If visible organic matter is present, use INTERCEPT at 2/3 oz./gal of water (0.5% concentration) **which is two full strokes of the hand-pump (2 oz.) to three gallons of water.** For best manual cleaning results, mix INTERCEPT with cool to warm water (20°C-35°C) (38°F – 95°F) and ensure a minimum contact time of (1) minute. Rinse all surfaces and internal channels thoroughly with water. For use in automated washers, follow washer manufacturer's recommendations. Do not mix INTERCEPT with bleach.

**PRECAUTIONARY STATEMENTS:** Use of personal protective equipment including appropriate clothing, gloves, and safety glasses must be worn. Refer to the material Safety Data Sheet (MSDS) for additional information.

**STATEMENT OF PRACTICAL TREATMENT:** In case of contact with eyes, flush with copious amounts of water for 15 minutes. In case of ingestion, DO NOT induce vomiting. Contact CHEMTREC for emergency medical treatment at 1-800-424-9300.

**STORAGE & DISPOSAL:** Store at room temperature 15°C-25°C (59°F-77°F). INTERCEPT has a two (2) year shelf life from date of manufacture. Discard unused solution according to local, state and federal regulations. DO NOT reuse container.

Read MSDS for additional precautionary statements

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© 2013 MEDIVATORS Inc. 67199-666/C
MEDIVATORS Inc.
Minneapolis, MN 55447 U.S.A.