# MEDIVATORS™ Bite Blocks

**REF** CATALOG NUMBER  
Standard Adult Bite Block #100411  
Pediatric Bite Block #100412

## INTENDED USE:

The bite block is to be used to prevent a patient from biting on the insertion tube of a flexible endoscope during upper GI endoscopic examinations.

## WARNINGS AND PRECAUTIONS:

- **DO NOT** attempt to use the bite block for any purpose other than that for which it has been designed.
- **NEVER** attempt to insert any endoscopic instruments or accessories with outside diameters exceeding 20mm through the adult bite block or 14.30mm through the pediatric bite block.
- **ALWAYS** detach both ends of the flexible strap from the tabbed ends of the bite block before attempting to remove the bite block from the patient's mouth.
- Carefully inspect all inner and outer surfaces of the product for cracked surfaces or any sharp edges. Do not use if product is damaged.
- Contact customer service for assistance.
- Federal (USA) law restricts the sale of this device by or on the order of a physician.

**NON STERILE**

## INSTRUCTIONS FOR USE:

### 100411 Standard Adult Bite Block and 100412 Pediatric Bite Block

1. Open the pouch and remove the bite block and attached flexible strap.
2. Place the bite block into the patient's mouth, in between the patient's teeth. Connect one end of the flexible strap through one of the side ventilation holes on the bite block and secure strap to bite block tab.
3. Position the unconnected end of the strap around the patient's head and insert it through the bite block side ventilation hole. Secure the strap to the remaining bite block tab. Adjust the strap to fit snugly, but do not over tighten.
4. When removing the bite block, always disconnect the flexible strap from both tabs of the bite block before attempting to remove the bite block from the patient's mouth. The bite block is for single use only, and should be disposed of properly after use.

Manufactured in the USA by:

Medivators Inc.  
3150 Pollock Drive  
Conroe, TX 77303  
Tel: +1.936.539.0391

Medivators Inc.  
14605 28th Avenue North  
Minneapolis, MN 55447 USA  
Toll Free: +1.800.444.4729

Medivators BV   
Sourethweg 11  
6422PC Heerlen  
The Netherlands  
Tel: +31.45.5.471.471

Cantel Medical Asia/Pacific Pte. Ltd.  
1A International Business Park  
#05-01 Singapore 609933  
Tel: +65.6227.9696

[www.medivators.com](http://www.medivators.com)

**NON STERILE**

SINGLE USE ONLY  
DO NOT REUSE

MEDIVATORS™ is a trademark of Medivators Inc.

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 **MEDIVATORS**  
A Cantel Medical Company

Manufactured in the USA by:

Medivators Inc.  
3150 Pollock Drive  
Conroe, TX 77303  
Tel: +1.936.539.0391

[www.medivators.com](http://www.medivators.com)

Medivators Inc.  
14605 28th Avenue North  
Minneapolis, MN 55447 USA  
Toll Free: +1.800.444.4729

Medivators BV   
Sourethweg 11  
6422PC Heerlen  
The Netherlands  
Tel: +31.45.5.471.471

Cantel Medical Asia/Pacific Pte. Ltd.  
1A International Business Park  
#05-01 Singapore 609933  
Tel: +65.6227.9696

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