Using SCOPE BUDDY™ PLUS Endoscope Flushing Aid 24-Hour Multi-Use Tubing Universal Hookup (SB-110-HU1016) with FUJIFILM™ 700 Series Endoscopes



## **Endoscope Connection**

FUJIFILM™ 700 Series Endoscopes require a unique connection method for compatibility with SCOPE BUDDY™ PLUS Endoscope Flushing Aid. The 24-Hour Multi-Use Tubing Hookup for FUJIFILM Endoscopes (SB-110-HU1009) should NOT be used to connect FUJIFILM 700 Series Endoscopes with the configurations shown in the Instruction for Use.

It is recommended to use the 24-Hour Multi-Use Tubing Universal Hookup (SB-110-HU1016) to connect to the FUJIFILM 700 Series Endoscopes. The Universal Hookup can be connected to FUJIFILM 700 Series Endoscopes using the FUJIFILM WA-007 and WA-010 cleaning adapters.

Please review the instructions and pictures below for information about this connection method.

1. Connect the WA-010 cleaning adapter to the control body portion of the endoscope and place the biopsy cap over the biopsy port.







2. Connect the WA-007 cleaning adapter to the light guide end of the endoscope.





3. Connect the red luer of the Universal Hookup (SB-110-HU1016) to the clear tube luer (suction channel) of the WA-007 cleaning adapter. Connect one of the clear luers of the Universal Hookup (SB-110-HU1016) to the blue tubing luer (air/water channel) of the WA-007 cleaning adapter





4. If the endoscope has a forward water jet, connect one of the clear luers to that channel. Clamp off any remaining unused tubes. The Universal Hookup is now ready to be connected to the SCOPE BUDDY™ PLUS Endoscope Flushing Aid pump tube.



If you have any questions or concerns regarding this product bulletin, please contact Customer Care at 1-800-444-4729.

SCOPE BUDDYTM is a trademark of Medivators Inc. FUJIFILM TM is a trademark of FUJILM Corporation.



Global Headquarters | p: +1.973.890.7220 Medical Division Headquarters | p: +1.763.553.3300

