

For advanced enzymatic cleaning, get MediClean EZ™



Endoscope reprocessing the way it should be™

MEDIVATORS®

Multi-enzymatic detergent with advanced proteolytic cleaning action

MediClean EZ™ Detergent...

- Provides fast enzymatic removal of blood, proteins, fats, carbohydrates, mucopolysaccharides and lipids
- Meets endoscope manufacturers' criteria for endoscope detergents
- Low foaming, neutral pH, medical grade detergent
- Cost-effective, concentrated formula reduces per-use cost
- Green tint for easy identification of dispensed detergent



Easy to use with a wide range of applications...

- Formulated for cleaning surgical instruments, immersible endoscopes, medical devices and accessories
- Supplied ready to use in a convenient dispensing container
- No premeasuring: Simply use 1 pump per 2 gallons of water
- For use in manual cleaning, automated washers and ultrasonic cleaners
- Environmentally responsible and biodegradable

MediClean EZ Ordering Information

Part number	Description
ML02-0131	Enzymatic Detergent (4 one-gallon bottles per case)



Minntech Corporation
14605 28th Avenue North
Minneapolis MN 55447
Tel: +1-763-553-3300
Fax: +1-763-553-3387

Minntech Asia/Pacific Pte Ltd
#05-01 Singapore 609933
Tel: +65-6227-9698
Fax: +65-6225-6848

Minntech BV
Sourethweg 11
6422 PC Heerlen
The Netherlands
Tel: +31-45-5-471-471
Fax: +31-45-5-429-695

Minntech China, Rep. Office
Room 708, 7th Floor Kaiheng Center,
No. 2 Chaoyangmennei Street,
Dongcheng District, Beijing China, 100010
Tel: +8610-6567-8446
Fax: +8610-6567-8445

For additional information please contact your Medivators representative
or visit our website at www.minntech.com/medivators