

Product Bulletin

For Distribution To:

- International Distributors
- International Management
- U.S. Field Sales & Clinical Specialists
- U.S. Technical & Field Service

Date: February 20, 2014

RE: DEFENDO™ SINGLE USE AIR WATER AND SUCTION VALVES

Dear Valued Customer:

The DEFENDO Single Use Air Water Valves, Suction Valves, Biopsy Valves, and Auxiliary Water connector for OLYMPUS® GI endoscopes (included in kits 100305, 100306, and 100310) have been internally qualified by MEDIVATORS and are considered compatible disposable replacement accessories for OLYMPUS GI endoscopes (140/160/180/190 Series).

DEFENDO Single Use OLYMPUS Accessories are packaged sterile, used for a single use, and then discarded. When used according to the Instructions for Use, the DEFENDO Single Use accessories provide quality alternatives to OLYMPUS reusable accessories (MH-438 Air Water valve, MH-443 Suction valve, and MB-358 Biopsy valve) while preventing potential cross contamination due to elimination of reprocessing steps.

Additionally, DEFENDO Single Use OLYMPUS Accessories do not cause any damage or harm to the scope ports. As the components used to create DEFENDO accessories are soft rubber and plastic, they cannot harm the metal ports into which they are inserted.

For further information, please refer to the DEFENDO accessory Instructions for Use.

A copy of this Product Bulletin is available on the MEDIVATORS website located at www.medivators.com. If you have any questions regarding this Product bulletin, please contact MEDIVATORS Customer Support at 1-800-444-4729.

DEFENDO™ is a trademark of MEDIVATORS Inc.
OLYMPUS® is a registered trademark of Olympus Corporation.



MEDIVATORS
3150 Pollok Dr.
Conroe, TX
77303 USA
Tel: +1.936.539.0391
Fax: +1.936.539.2381

MEDIVATORS Inc. Headquarters
14605 28th Avenue North
Minneapolis, MN
55447-4822 USA
Tel: +1.800.444.4729
Fax: +763.553.3387

MEDIVATORS BV
Sourethweg 11
6422PC Heerlen
The Netherlands
Tel: +31.45.5.471.471
Fax: +31.45.5.429.695

www.medivators.com