

# Product Bulletin

## For Distribution To:

- International Distributors
- International Management
- U.S. Field Sales & Clinical Specialists
- U.S. Technical & Field Service

**Date: March 27, 2014**

## **RE: DEFENDO™ SINGLE USE SUCTION VALVE FOR OLYMPUS® GI ENDOSCOPES**

Please make note of the product improvement outlined below:

Application of a thin layer of lubrication to the disposable DEFENDO™ Single Use Suction Valve stem to aide in ease of valve actuation.

This improvement is applicable to the following product lines and will be first implemented in the lot numbers indicated:

100305 – Lot # 738450 Rev P.

100306 – Lot #739882 Rev. M

100310 – Lot #739030 Rev. E

\*All 800 series procedure kits and custom procedure kits will include the following:

C305 – Lot # 739029 Rev. G

C310 – Lot # 739032 Rev. F

Please contact Erica Dodson Imhoff with any questions or concerns – [edodson@medivators.com](mailto:edodson@medivators.com) or 936-521-3983.

DEFENDO™ is a trademark of Medivators Inc.  
OLYMPUS® is a registered trademark of Olympus Corporation.



MEDIVATORS  
3150 Pollok Dr.  
Conroe, TX  
77303 USA  
Tel: +1.936.539.0391  
Fax: +1.936.539.2381

MEDIVATORS Inc. Headquarters  
14605 28th Avenue North  
Minneapolis, MN  
55447-4822 USA  
Tel: +1.800.444.4729  
Fax: +763.553.3387

MEDIVATORS BV  
Sourthweg 11  
6422PC Heerlen  
The Netherlands  
Tel: +31.45.5.471.471  
Fax: +31.45.5.429.695

[www.medivators.com](http://www.medivators.com)