


# PULL THRU™

Single Use Channel Cleaning Device



≈ 2200 mm

2.8 - 5.0 mm id



Single use only

Use wet

Read Instructions before use

see [www.pullthru.net.nz](http://www.pullthru.net.nz) for more information





MEDIVATORS  
3150 Pollock Drive  
Conroe, TX 77303  
Tel: +1.936.539.0391  
[www.medivators.com](http://www.medivators.com)

MEDIVATORS Inc. HQ  
14605 28th Avenue North  
Minneapolis, MN 55447 USA  
Toll Free: +1.800.444.4729

MEDIVATORS BV  
Sourethweg 11  
6422PC Heerlen  
The Netherlands  
Tel: +31.45.5.471.471

EC REP

New Zealand Patent Number 540125, Other International Patents Pending  
Copyright © 2012 F R Galantai Manufacturing Limited  
Trademarks are the property of their respective owners




Distal end

Y = 1.4-2.6 mm id  
**PULL THRU™**


IFU-67 /REV A

AA




A1

Pentax®




A1

Fujinon®




A2

Pentax®



A2


Fujinon®



A3


B1

Pentax®




B1


Olympus®



C1



B2




C2


A3

B3

C3



Pentax® Endoscopes 30, 30k/40 series



Pentax® Endoscopes 30, 30k/40 series

