


ENDO SMARTCAP[®] ENDOGATOR[®]

Disposable Tubing | **Procedure**





THE COMPLETE CIRCLE OF PROTECTION

As the global vanguard in infection prevention, only Cantel Medical delivers the Complete Circle of Protection, a full-value, proactive partnership dedicated to helping you remove risk, streamline operational efficiencies and optimize your success.

Disposable GI Endoscopy Irrigation Solutions designed with patient safety in mind.

Infection Prevention

- No need for reprocessing
- Reduced risk of cross contamination

Efficiency

- Minimizes reprocessing requirements
- Creates consistent practices

Compatibility

- Use of CO2 or air insufflation
- Works with all irrigation pumps and sterile water bottles





PROCEDURE

Reducing the risk of patient cross contamination is at the forefront of infection prevention. Cantel innovates infection control products designed to improve patient outcomes, while increasing procedural efficiency.

ENDO SMARTCAP® ENDO GATOR®

Disposable Irrigation Tubing Solutions | **Procedure**



ENDO SMARTCAP® Disposable Tubing

- Disposable alternative to the reusable water bottle
- Increases infection prevention measures



ENDO GATOR® Disposable Irrigation Tubing

- Cost effective and completely disposable
- 24-hour use without the need to reprocess
- Increases infection prevention measures



ENDO GATOR™ Disposable Irrigation Tubing Kit

- One order number
- Innovative lens cleaning/irrigation system
- Option to use room air or CO₂ insufflation

PENTAX®
Single Use
Connector



OLYMPUS®
Single Use
Connector



ENDO GATOR® Disposable Irrigation Hybrid Tubing

- Combined lens cleaning and irrigation system
- Uses room air or CO₂
- Fits any sterile water bottle

ENDO SMARTCAP® Disposable Tubing

Ordering Information

Model Numbers		Description	Units/Box
US	OUS		
100145	100145U	Disposable Tubing – For OLYMPUS® 140/160/180 Series GI Endoscopes (CO ₂ Also Available)	10
100150	100150U	Disposable Tubing – For OLYMPUS 100/130 Series GI Endoscopes (CO ₂ Also Available)	10
100160	100160U	Disposable Tubing – For PENTAX® GI Endoscopes (CO ₂ Also Available)	10
100165	100165U	Disposable Tubing – For FUJIFILM™ GI Endoscopes (CO ₂ Also Available)	10
100141	100141	Adapter – For FUJIFILM 530 Series Colonoscopes	50

ENDOGATOR® Disposable Irrigation Tubing

Ordering Information

Model Numbers		Description	Units/Box
US	OUS		
100130	100130	Disposable Irrigation Tubing – For ENDO STRATUS®, OLYMPUS OFF, or ENDOSTAT™ II Irrigation Pumps	10
200230	200230U	Disposable Irrigation Tubing – For OLYMPUS AFU-100, OLYMPUS OFF, ERBE® EIP2 or EGP-100 Irrigation Pumps	10
200230F	200230FU	Disposable Irrigation Tubing – For EGP-100 Irrigation Pump with FUJIFILM Connector	10
100241	100241	Single Use Auxiliary Water Jet Connector – For OLYMPUS 140/160/180 series GI Endoscopes	50
100242	100242	Single Use Auxiliary Water Jet Connector – For PENTAX GI Endoscopes	50
100130F	100130FU	For ENDO STRATUS, OLYMPUS OFF, or ENDOSTAT II Irrigation Pumps for FUJIFILM GI Endoscopes	10
100135	100135	Channel Irrigation Tubing	10

ENDOGATOR® Disposable Hybrid Tubing

Ordering Information

Model Numbers		Description	Units/Box
US	OUS		
100605	100605S	Hybrid Tubing with CO ₂ connection – For OLYMPUS 140/160/180/190 Series GI Endoscopes with ENDO STRATUS, OLYMPUS OFF or ENDOSTAT II Irrigation Pumps with 24 hour auxiliary water jet connector	10
100606	100606S	Hybrid Tubing with CO ₂ connection – For OLYMPUS 140/160/180/190 Series GI Endoscopes with OLYMPUS AFU-100, OLYMPUS OFF, EGP-100 or ERBE® EIP2 Irrigation Pumps with 24 hour auxiliary water jet connector	10
100608	100608S	Hybrid Tubing with CO ₂ connection – For PENTAX GI Endoscopes with EGP-100 or ERBE EIP2 Irrigation Pumps	10
100609	100609S	Hybrid Tubing with CO ₂ connection – For OLYMPUS 140/160/180/190 Series GI Endoscopes with ENDO STRATUS, OLYMPUS OFF or ENDOSTAT II Irrigation Pumps	10
100610	100610S	Hybrid Tubing with CO ₂ connection – For OLYMPUS 140/160/180/190 Series GI Endoscopes with OLYMPUS AFU-100, OLYMPUS OFF	10

ENDOGATOR™ Tubing Kit Ordering Information

Please contact your MEDIVATORS sales consultant for kit option

ENDO SMARTCAP®, ENDOGATOR® AND ENDO STRATUS® are registered trademarks of Medivators Inc.
 PENTAX® is a registered trademark of Hoya Corporation. FUJIFILM™ is a trademark of FUJIFILM Corporation.
 OLYMPUS® is a registered trademark of Olympus Corporation. ERBE® is a registered trademark of ERBE Elektromedizin GmbH.
 ENDOSTAT™ is a trademark of Boston Scientific Corporation.



**TO PLACE
AN ORDER**

Call: 1.800.328.3340 (Customer Service)
 Fax: 1.800.686.8493
 Email: custserv@medivators.com

www.medivators.com

© 2017 Medivators Inc. P/N 50098-246D