

INTERCEPT®

Detergent Family | **Manual Clean**





THE COMPLETE CIRCLE OF PROTECTION

As the global vanguard in infection prevention, only Cantel Medical delivers the Complete Circle of Protection, a full-value, proactive partnership dedicated to helping you remove risk, streamline operational efficiencies and optimize your success.

Solutions for Superior Cleaning

INTERCEPT® Detergent

Effective reprocessing procedures start with manual cleaning. INTERCEPT Detergent is a unique surfactant based, non-enzymatic formulation specifically developed for manual or automated cleaning of endoscopes, surgical instruments, and accessories. The concentrated solution has a fast contact time and low use concentration to save your facility time and money. INTERCEPT Detergent meets society guidelines for effective manual cleaning.

- Clinically proven to provide superior removal of biological and organic soils found on instruments
- Effectively aids in the removal of biofilm
- Low foaming and neutral pH formulation is tough on organic soils, while safe on instruments
- Cost-effective, concentrated formula reduces per-use costs
- Fast one-minute contact time increases reprocessing efficiencies
- Environmentally responsible and biodegradable chemistry

INTERCEPT® Wipes

INTERCEPT Wipes provide superior cleaning of instrument surfaces at the bedside or during manual cleaning. The ready-to-use wipes are pre-saturated with clinically tested INTERCEPT Detergent for convenient and cost effective removal of soils from endoscopes, instruments and counter tops.

- Lint-free, durable and highly absorbent wipes provide superior cleaning of instrument surfaces
- Large 7"x10" wipe increases surface area and delivers a ready-to-use detergent solution in every sheet
- The low cost-per-use standardizes practices and saves money.
- Eliminates the need for costly sponges during manual cleaning



INTERCEPT Wipes

INTERCEPT Detergent



MANUAL CLEAN

Manual cleaning is essential in removing bioburden and preparing the endoscope for high-level disinfection. For optimal results, Cantel's manual cleaning products and detergents are designed to be a critical first step to efficiently and effectively eliminating bacteria on endoscopes prior to reprocessing.

INTERCEPT®

Detergent Family | **Manual Clean**

Endo Procedure Kits

The essential you need for safe and effective endoscopy procedures.
(Multiple varieties available; consult your Sales Representative for details)



INTERCEPT™ Bedside Kit

Arm your staff with INTERCEPT Bedside Kits for a convenient and superior precleaning process.

- Pre-measured INTERCEPT Detergent packets ensure proper dilution and improves the removal of organic debris at bedside
- Universal basin offers flexibility for 250 mls or 500 mls
- Sturdy basin design reduces spill potential
- Convenient, disposable solution eliminates contamination risks associated with reusing towels and basins



Easy-to-use sponge locks in fluid for enhanced wiping of the insertion tube

Using INTERCEPT® Detergent Throughout Your Facility will Increase Compliance with Manual Cleaning Standards

ASGE

“To date, all published occurrences of pathogen transmission related to GI endoscopy have been associated with failure to follow established cleaning and disinfection/sterilization guidelines or use of defective equipment.”

2011 ASGE Multisociety Guideline on Reprocessing Flexible Gastrointestinal Endoscopes

CDC

“Thorough cleaning is required before high-level disinfection and sterilization because inorganic and organic materials that remain on the surfaces of instruments interfere with the effectiveness of these processes. Also, if soiled materials dry or bake onto the instruments, the removal process becomes more difficult and the disinfection or sterilization process less effective or ineffective. Surgical instruments should be presoaked or rinsed to prevent drying of blood and to soften or remove blood from the instruments.”

2008 CDC Guidelines for Disinfection and Sterilization in Healthcare Facilities

SGNA

“Complex endoscope design features may allow organic debris and microorganisms to accumulate, making manual cleaning essential. Biofilm formation may harbor microorganisms, making strict and meticulous adherence to reprocessing guidelines imperative in order to prevent cross-contamination between patients or hospital-acquired infections. Prompt efficient cleaning processes are the best defense against biofilm formation”

2012 SGNA Standards of Infection Control in Reprocessing of Flexible Gastrointestinal Endoscopes

INTERCEPT Detergent Family Ordering Information

Part Number	Description
ML02-0106	INTERCEPT Detergent (4 one-gallon bottles/case)
ML02-0107	INTERCEPT™ Wipes (4 containers/case, 100 wipes/container)
100900	INTERCEPT™ Bedside Kit (box of 50 kits containing detergent packet, sponge and disposable basin)
Consult Your Sales Representative	Endo Procedure Kits



INTERCEPT® and DEFENDO® are registered trademarks of Medivators Inc.

CANTEL
Global Headquarters
150 Clove Road - 9th Floor
Little Falls, NJ 07424-2139
Phone: +1.973.890-7220

MEDIVATORS
Global Headquarters
14605 28th Avenue North
Minneapolis, MN 55447-4822
Phone: +1.763.553.3300 Fax: +1.763.553.3387

TO PLACE
AN ORDER

Call: 1.800.328.3340 (Customer Service)
Fax: 1.800.686.8493
Email: custserv@medivators.com

www.medivators.com

© 2017 Medivators Inc. P/N 50098-669C