



# Certificate of Analysis

## RAPICIDE® PA

### Part B

Product Code: ML02-0148  
Lot Number: 00017619  
Expiration Date: 2020/10

Components	Specification	Analysis
Sodium	743.0 - 909.0 mEq/L	775.4 mEq/L
Benzotriazole	0.9 – 1.1% w/w	1.0 %w/w
pH	11.4 to 12.4 as is	11.7

This lot has been reviewed and meets specified requirements.

Quality Analyst

Analysis Date

30/04/19

Quality Assurance

Release Date

02-MAR-19

**Medivators Inc.**  
14605 28th Avenue North  
Plymouth, Minnesota 55447  
Telephone: (763) 553-3300  
Fax (763) 553-3387

RAPICIDE® is a registered trademark of Medivators Inc.