




MEDIVATORS
A Cantel Medical Company


Certificate of Analysis
RAPICIDE® PA
Part A

Product Code: ML02-0143
 Lot Number: 00020217
 Expiration Date: 2019/01

Components	Specifications	Analysis
H ₂ O ₂	20.0 – 24.0% (w/w)	23.5 %
PAA	4.0 – 6.0 % (w/w)	5.6 %
Acetic Acid	8.0 – 10.0% (w/w)	8.4 %
pH, (as is)	0.5 – 1.1 (as is)	0.5

This lot has been reviewed and meets specified requirements.


 _____ 04 July 2017
 Quality Analyst Analysis Date


 _____ 04 July 2017
 Quality Assurance Release Date

Medivators Inc.
 14605 28th Avenue North
 Minneapolis, MN 55447 USA
 Telephone: (763) 553-3300
 Fax (763) 553-3387
 www.medivators.com

RAPICIDE® is a registered trademark of Medivators Inc.