

# Certificate of Analysis


## RAPICIDE<sup>®</sup> PA

### Part A


Product Code: ML02-0147  
 Lot Number: 00023019  
 Expiration Date: 2020/12

Components	Specifications	Analysis
H <sub>2</sub> O <sub>2</sub>	20.0 – 24.0% (w/w)	23.5 %
PAA	4.0 – 6.0 % (w/w)	5.0 %
Acetic Acid	8.0 – 10.0% (w/w)	9.2 %
pH, (as is)	0.5 – 1.1 (as is)	1.0

This lot has been reviewed and meets specified requirements.

  
 \_\_\_\_\_  
 Quality Analyst

10/26/19  
 \_\_\_\_\_  
 Analysis Date

  
 \_\_\_\_\_  
 Quality Assurance

13-JUN-19  
 \_\_\_\_\_  
 Release Date

**Medivators Inc.**  
 14605 28th Avenue North  
 Minneapolis, MN 55447 USA  
 Telephone: (763) 553-3300  
 Fax (763) 553-3387  
 www.medivators.com

RAPICIDE<sup>®</sup> is a registered trademark of Medivators Inc.