



# Certificate of Analysis

## INTERCEPT<sup>®</sup> PLUS

### Detergent & Cleaner


Part No. ML02-0145

Lot Number: 00032717


Expiration Date: 2019/10

Property	Test Method	Specification	Analysis
Appearance	Visual	Clear Non-viscous Liquid	Pass
pH (release) at 25°C	20002-154	11.5 – 12.5	12.0
Phosphate Level	20602-105	80.83 – 98.79 g/L	92.82 g/L

This lot has been reviewed and meets specified requirements.

  
\_\_\_\_\_  
Quality Analyst

17 October 2017  
Analysis Date

  
\_\_\_\_\_  
Quality Assurance

23-oct-17  
Release Date

Medivators Inc.  
14605 28th Avenue North  
Minneapolis, MN 55447 USA  
Telephone: (763) 553-3300  
Fax (763) 553-3387  
www.medivators.com

INTERCEPT<sup>®</sup> is a registered trademark of Medivators Inc.