

Certificate of Analysis
RAPICIDE[®] PA
Part A

Product Code: ML02-0147
Lot Number: 00035718
Expiration Date: 2020/04

Components	Specifications	Analysis
H ₂ O ₂	20.0 – 24.0% (w/w)	<u>23.4 %</u>
PAA	4.0 – 6.0 % (w/w)	<u>5.9 %</u>
Acetic Acid	8.0 – 10.0% (w/w)	<u>8.5 %</u>
pH, (as is)	0.5 – 1.1 (as is)	<u>1.0</u>

This lot has been reviewed and meets specified requirements.

FDC

Quality Analyst

31-OCT-18

Analysis Date

[Signature]

Quality Assurance

31-OCT-18

Release Date

Medivators Inc.
14605 28th Avenue North
Minneapolis, MN 55447 USA
Telephone: (763) 553-3300
Fax (763) 553-3387
www.medivators.com

RAPICIDE[®] is a registered trademark of Medivators Inc.