

Certificate of Analysis

RAPICIDE® PA

Part A


Product Code: **ML02-0143**

Lot Number: 00039017

Expiration Date: 2019/05

Components	Specifications	Analysis
H ₂ O ₂	20.0 – 24.0% (w/w)	<u>22.1 %</u>
PAA	4.0 – 6.0 % (w/w)	<u>5.6 %</u>
Acetic Acid	8.0 – 10.0% (w/w)	<u>8.2 %</u>
pH, (as is)	0.5 – 1.1 (as is)	<u>0.5</u>

This lot has been reviewed and meets specified requirements.



Quality Analyst

30/11/2017

Analysis Date



Quality Assurance

30/11/17

Release Date

Medivators Inc.
14605 28th Avenue North
Minneapolis, MN 55447 USA
Telephone: (763) 553-3300
Fax (763) 553-3387
www.medivators.com

RAPICIDE® is a registered trademark of Medivators Inc.